MyPaclitaxel™ Calibrator Kit


This package insert must be read carefully prior to product use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Customer Service:

Telephone: +1-610-419-6731

Fax: +1-484-547-0590

Key to Symbols Used

| | | | |
|---|--|--|--|
|  | <i>in vitro</i> Diagnostic Device |  | Temperature Limitation |
|  | Batch Code (Lot) |  | Use By |
|  | Catalog Number |  | Authorized Representative in the European Community |
|  | Caution: Consult accompanying documents |  | Manufacturer |
|  | Consult Instructions for Use |  | Calibrator A through Calibrator F |

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

Saladax Biomedical, Inc.
116 Research Drive
Bethlehem, PA 18015 USA
www.saladax.com

Intended Use

The Saladax MyPaclitaxel Calibrator Kit is intended for use in calibration of the Saladax MyPaclitaxel Assay.

Contents

The MyPaclitaxel Calibrator Kit consists of a buffer matrix with the following concentrations of paclitaxel:

| Vial | | Concentration (ng/mL) | Quantity | Volume |
|------|---|--------------------------|----------|--------|
| CAL | A | 0 | 1 | 3 mL |
| CAL | B | 40 | 1 | 3 mL |
| CAL | C | 80 | 1 | 3 mL |
| CAL | D | 160 | 1 | 3 mL |
| CAL | E | 240 | 1 | 3 mL |
| CAL | F | 320 | 1 | 3 mL |

Standardization

There is no internationally recognized standard for paclitaxel. The MyPaclitaxel Calibrators are prepared gravimetrically by dilution of USP Paclitaxel into a buffer matrix at concentrations of 0, 40, 80, 160, 240, and 320 ng/mL.

Saladax Biomedical, Inc. maintains master calibrators that have been confirmed by a reference method. MyPaclitaxel Calibrators are traceable to this master set of calibrators.

Precautions and Warnings

For In Vitro Diagnostic Use.

The calibrators in this set are designed for use as a unit. Do not substitute or mix calibrators with those from other lots.

Exercise the normal precautions required for handling all laboratory reagents.

All MyPaclitaxel Calibrators contain less than 0.1% sodium azide. Avoid contact with skin and mucous membranes. Flush affected areas with copious amounts of water. Seek immediate medical attention if calibrators are ingested or come into contact with eyes.

Storage and Stability

Store calibrators refrigerated at 2-8°C (35-46°F). Do not freeze. Calibrators are stable until the expiration date when stored as directed. Once opened calibrators are stable for 2 months when handled as directed.

Limitations of the Procedure

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, calibrators, correct product storage, and good laboratory technique.

Instructions for Use

- Refer to the Saladax MyPaclitaxel Assay package insert, included in the reagent kit, for a complete summary and explanation of the test.
- Use calibrators immediately upon removal from 2-8°C storage.
- Mix each calibrator by gentle inversion several times before dispensing.
- After each use, tightly close the caps and return calibrators to 2-8°C.

Bubbles may interfere with proper detection of calibrator level in the sample cup, causing insufficient calibrator aspiration that could impact results.

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