



MyCare™ INSITE

SYMBOLS

	Consult Instructions for Use		Electronic Equipment: Do not dispose with household waste
	in vitro Diagnostic Device		Manufacturer
	Temperature Limitation		Serial number
	Catalog Number		Caution
	Electrical Symbol		Authorized Representative in the European Community
	Near Patient Testing		Indicates conformity with the provisions of the applicable EU directives and/or regulations
	Date of Manufacture		Relative humidity limitation

INTENDED PURPOSE

The MyCare Insite is intended as an automated general-purpose laboratory instrument for the quantitative determination of clinical wet chemistry parameters. The Analyser is intended to be used solely in conjunction with MyCare Insite Test kits and is designed to be used either in a clinical laboratory or in a near-patient setting by health-care professionals.

PRINCIPLE OF OPERATION

The MyCare Insite is used to quantify an analyte concentration in samples with an optical absorbance measurement. To perform a test the MyCare Insite Analyser is loaded with a test cartridge containing reagents and the sample to be tested.

The Analyser is equipped with an RFID card-reader module. RFID cards are necessary for performing any testing procedures. They are included in the test kits and contain all the specific steps for the various tests, the lot data, as well as the calibration data. The Analyser performs the tests automatically according to the information on the RFID card. The test measured is dictated by the RFID card placed on the Analyser. RFID cards and tests are color-coded. Numerous types of tests can be performed automatically.

The sample and the reagents are automatically mixed within the Analyser. The photometer unit performs the analysis with a light emitting diode (LED). The absorption of light rays is determined during this process and the measured values are then automatically converted into the test result using mathematical methods.

The test cartridge design is self-contained so that the Analyser has no contact with reagent or sample.

The result is displayed on the Touch Screen. Optionally, results can be exported to an external computer or an HIS/LIS and can also be printed to an external printer.


After the test process the door opens automatically and the test cartridge can be removed and discarded. After confirming the result on the Touch Screen and closing the door, the Insite is ready to perform the next analysis.

With LED technology the Insite requires no maintenance.

ANALYSER COMPONENTS

The MyCare Insite consists of the MyCare Insite Analyser laboratory photometer and the MyCare Insite Touch Screen. The Insite package contains:

- User Manual
- MyCare Insite Analyser
- Power adapter
- Power cable
- Test Rack
- USB cable
- Touch Screen
- Touch Screen Stand
- Touch Screen Mounting instructions

Before first use remove the transport lock from the door of the Analyser.

Note: Keep both the transport lock and the original packaging. In case of a defect, the device must be returned in the original packaging with the transport lock installed to prevent damage in transport.

The Insite Technical Guide and videos for unpacking and setting up the Insite are at MyCareInsite.com. With WiFi Connectivity tap  on the Main Menu to access the information.

WARNINGS AND PRECAUTIONS

- For in vitro Diagnostic Use Only
- The Analyser performs self checks during measurements. Any electro-mechanical error will stop a result from being generated. If an error message appears on the screen, please consult the “ERROR MESSAGES AND POSSIBLE CAUSES” on page 10. Repeat the test.
- The door prevents ambient light, dust, dirt, and humidity from entering the Analyser during the testing process.
- Remove the test cartridge from the Analyser after every test. Dispose of the test cartridge appropriately.
- Make sure the door is closed whenever the Analyser is not in use.
- Do not attempt to open the Analyser door manually.
- Use only your fingertips to operate the Touch Screen! Do not use pens or other objects that may scratch or damage the screen.
- Do not use liquid cleaners on the Analyser or Touch Screen. Use a lint-free, slightly damp cloth only.
- Connect optional equipment only when the Analyser is unplugged. Please note that attaching optional equipment (e.g., a printer) can increase the amount of leakage current.
- Use ONLY manufacturer-approved tests. Otherwise, serious damage to the Analyser and/or inaccurate test results may occur.
- To ensure samples are correctly identified run only one sample at a time.
- To ensure safe operation, always use the MyCare Insite according to the instructions in the package insert.
- Avoid electromagnetic radiation and electrostatic discharge.
- Prevent physical access by unauthorized users
- Enable screen lock setting “Password” on the Touch Screen.
- Any serious incident that occurs in relation to the MyCare Insite shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- Residual Risks: User error may cause display of a wrong measurement result, delayed result, no result; and unauthorized access to information.

DESCRIPTION OF THE MYCARE INSITE

This section introduces the general characteristics of the MyCare Insite Analyser and Touch Screen.


RFID Card Well – The RFID card is placed here.

Indicator Lights – There are 3 LEDs on the front of the Analyser.

Door – When the door opens the test cartridge (Cuvette with Reagent Cap) can be inserted or removed from the inner Chamber.


Note: *do not attempt to open the door manually.*

Touch Screen – User interface to operate the Analyser & process the results.

MyCare Insite Analyser Interfaces

MyCare Analyser Back

Indicator Lights

MyCare Insite Analyser Front

LEDs on the MyCare Insite Analyser

 **This green LED indicates the Analyser status**

Flashing: the Analyser is warming up and testing cannot be started yet. (Note: warm-up takes approximately 10 min.)

On (green): the Analyser is warmed up and ready for use

Off: the Analyser is not plugged in

 **This blue LED indicates the Touch Screen connection status**

On (blue): connection with the Touch Screen is established

Off: the Touch Screen is not connected

 **This blue LED indicates Bluetooth connection status**

On (blue): connection with the Touch Screen is established

Flashing: Analyser Bluetooth is activated but no Touch Screen is connected via Bluetooth

Off: Bluetooth is switched off

Note: Use of the USB cable connection is recommended to keep the Touch Screen powered.

Operating and Storage Conditions

Operating conditions:

- 20 - 28°C; relative humidity: 10 – 85%
- A dry, clean, level surface with at least 10 cm of clearance on each side

Storage/Transport:

- 0 - 50°C; relative humidity: 5 – 85%
- Always use the transport lock when transporting the analyser or damage may result

SETTING UP THE MYCARE INSITE ANALYSER AND TOUCH SCREEN

AVOID

- Condensing humidity and water
- Heat and large temperature fluctuations
- Electromagnetic radiation
- Electrostatic discharge
- Vibrations (e.g. from centrifuges and dishwashers)
- Direct sunlight

Place the Analyser on a dry, clean, stable, and level surface.

Make sure the Analyser has at least 10 cm of table surface and clearance on each side, and that the Analyser can be easily disconnected from the power source. Allow the Analyser to acclimate to the ambient room temperature before operation.

Plug the Power Adapter (AC adapter plug) into the Analyser and the Power Cable into the wall socket.

The MyCare Insite is operated with a Touch Screen.

To ensure that the Analyser is properly grounded, use a 3-pole power outlet that includes neutral, hot, and ground sockets. The MyCare Insite must be connected to a power outlet with a 3-pole plug.

A certified electrician can confirm that the grounding meets requirements.

Ensure the Analyser is plugged in. Connect the Touch Screen to the USB port of the Analyser and then turn on the Touch Screen.

The MyCare Insite Launch Screen


The launch screen appears by pressing the Touch Screen home button or when the Touch Screen is not connected to the Analyser.

Touch the MyCare Insite logo to start the MyCare Insite application.

Touch Browser to open the Android web browser.

Touch Settings to open the Android settings menu.

1. Connection Screen

The connection screen appears when the Touch Screen is connected to the Analyser.

Tap Next to continue.

2. Main Menu

The MyCare Insite application starts and the main menu appears.

Note: while the MyCare Insite is warming-up (indicated by the flashing green LED on the front, see page 5) the "Measurement" button is grey and inactive.

How to Wake-Up the MyCare Insite or Switch It OFF

When the Touch Screen displays the main menu, the “power save” function dims the screen. Pressing the lock switch will re-illuminate the display to its customary level of brightness. It is not necessary to switch the MyCare Insite off every day. Turning off the device completely can only be done by disconnecting it from the power supply.

Materials Provided:

REF

MCI-EUR – MyCare Insite (consisting of MyCare Insite Analyser laboratory photometer and the MyCare Insite Touch Screen)

Materials Required – Provided Separately:

MyCare Insite Tests and Controls

Materials Not Provided – Optional:

- An external printer – Seiko Label Printer 650SE, serial 100 – 240 VAC
- An external barcode scanner - Honeywell Voyager 1200g, 3.8~13.2 cm (5 mils), 0.09 mm (3.5 mils), 100 – 240 VAC, for optional entry of patient or operator ID
- A PC – for the transfer of test data into a hospital or laboratory Information system (HIS/LIS)


ANALYSING A PATIENT SAMPLE OR CONTROL

Please refer to package insert provided for the MyCare Insite Test at MyCareInsite.com and follow all instructions including:

- The proper temperature of a test before it can be used
 - The exact volume of sample
 - The proper storage of the tests
-
- Be sure to handle the tests according to the instructions in the package inserts.
 - Be sure the test cartridge (Cuvette with Reagent Cap) is properly sealed before inserting into the Analyser.
 - Be sure the test cartridge (Cuvette with Reagent Cap) is fully inserted into the Analyser Chamber.
 - Use ONLY manufacturer-approved tests. Otherwise, serious damage to the Analyser and/or inaccurate test results may occur.
 - Do not attempt to open the door manually – open the door from the Service/User Service Menu.
 - A quality control program should be performed on a regular basis to verify the MyCare Insite is working properly and providing reliable results.
 - Data integrity can only be assured when controls and GCLP practices are used routinely.
 - The frequency of performing QC shall comply with local, state and/or national regulations.

1. RFID Card Placement

Take the provided RFID card out of the test kit and place it on the RFID card well on top of the Analyser.

3. Test Screen


The MyCare Insite displays the test name and number of tests remaining on the RFID card.

- Ensure the test name is the test being run
- Ensure the RFID and Test components have the same color coding

Touch    to proceed.

Touch  to return to the main menu.

5. Insert Cartridge and Close Door


The door opens automatically. Insert the cartridge (Cuvette with Reagent Cap) into the Analyser and close the door. The test starts automatically.

2. Main Menu

To start a measurement, select  in the Main Menu. The MyCare Insite automatically reads the RFID card and the Test Screen appears.

4. Input Menu

Select the desired input field (e.g. "Name"); the Touch Screen keyboard appears. Enter the information and confirm with "Done". Next, tap  and select an entry from the drop-down menu. Touch  to confirm the input or  to exit to the main menu.

6. End of Measurement


At the end of the measurement the door opens and the result appears on the Touch Screen. Remove the test cartridge and close the door for the next analysis.

VIEWING AND PROCESSING TEST RESULTS

After processing a test, the result screen appears on the Touch Screen:

1. Result Screen

The test result is displayed on the Touch Screen.

Use and to scroll back and forth within the results.

Tap to export to a HIS/LIS or PC or send as email.

Tap to print.

Tap to open the photometric data curve.

Tap to synchronize the result with MyCOS (MyCare Online Services) server.

Tap to delete the result.

Tap to exit to the main menu.

All test results are stored on the Touch Screen. They can be accessed as follows:

1. Main Menu

Tap in the main menu to open the Result List.

2. Result List Screen

The stored results can be browsed by swiping up and down.

Tap a result to open the result screen (see above), to return to the result list from a result screen press the Return Arrow .


Tap to open the filter options.

Tap to export or mail.

Tap to print, mail, or export.

Tap to synchronize results.

2.1 Filter Options

Results can be filtered by date, test, and Patient ID.

Tap to set the desired start date, and to set the desired end date.

Tap to select a test name in the drop-down menu (e.g., CLZ).

Tap to filter by patient ID.

ERROR MESSAGES AND POSSIBLE CAUSES


Door opening during a measurement indicates an error. Consult the following error messages. If an error occurs during the measurement, the test must be repeated.

Error message	Cause	Correction
Invalid card	A wrong, defective, or expired RFID card has been used. Incompatible application. A defective RFID module.	Use a new test kit RFID card. Update Insite application. Please contact your Official Saladax Distributor.
Tests expired	The test has passed its expiration date.	Use a new test that has not expired.
Wrong date set	The date set on the Touch Screen is wrong.	Change the date to "today" as described in the Touch Screen Starter Manual.
Door blocked	The test Cartridge is blocking the door because the Cartridge has not been inserted completely or the Reagent Cap does not close the Cuvette firmly enough.	Reposition the test cartridge or tighten the Reagent Cap on top of the Cuvette. Please contact your Official Saladax Distributor if the door fails to open at the end of a run.
Wrong cap Missing cap Missing cartridge	The wrong Reagent Cap is being used, or the Reagent Cap is missing, or the Cartridge is missing, or an already used Cartridge is inserted.	Use the correct Cuvette and Reagent Cap.
Bolt blocked	The test Cuvette blocks the bolt (that locks the Chamber) because the wrong Reagent Cap is being used.	Use the correct Cuvette and Reagent Cap.
Measurement overflow	The photometric measurement value lies outside the measuring range (e.g. a cold Cuvette and Reagent Cap have been used).	Repeat the test using a new Cuvette and Reagent Cap.
Blank error	The photometric measurement value is outside the measuring range (without the cartridge).	Repeat the test after restarting the Analyser.
Temperature error	The temperature is outside the range.	Repeat the test after restarting the Analyser.
Mix error	Analyser fails to perform test due to possible hardware error.	Please contact your Official Saladax Distributor.
Sample volume error	The first measurement (before the reaction start) is above the specified limit.	Repeat the test using a new Cuvette and Reagent Cap. Ensure the Cuvette and Reagent Cap have warmed up for at least 10 minutes and there is no condensation on the Cuvette. In the case of condensation wipe the Cuvette with a lint free cloth.
Cuvette temperature error	The Cuvette's temperature is too low for a proper test.	Repeat the test and refer to the package insert for proper Cuvette use.
Reaction error	The Cuvette may have been inserted without a Reagent Cap or there was another issue with the Reagent Cap.	Repeat the test using a new Cuvette and Reagent Cap.

CLEANING INSTRUCTIONS

Cleaning the Door & Chamber

- Clean the door with a clean, lint-free cloth moistened with isopropyl alcohol. Wipe down the surface.
- The cartridge area may be cleaned with a cotton swab moistened with isopropyl alcohol.

Clean the Touch Screen display with a clean, lint-free cloth moistened with water

- Do not use any liquid other than water when cleaning the Touch Screen!
- To avoid damage avoid excess water on the Touch Screen!
- Never use isopropyl alcohol on the Touch Screen!

Cleaning the Exterior

- The surface of the MyCare Insite Analyser should be cleaned with a clean, lint-free cloth only moistened with isopropyl alcohol. Wipe down the exposed surfaces.
- Follow site protocols for cleaning intervals.

TECHNICAL SPECIFICATIONS

MyCare Insite Analyser

Photometer resolution	0.0001 ABS
Reproducibility	< 1.5% CV at 1 OD
Linearity	0.1000 – 3.0000 OD better than $\pm 1.5\%$ and ± 0.01 OD
Temperature control	Electrical temperature control of the photometer unit to $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$
Fuse	2.5 A, self-healing
Dimensions	160 x 135 x 145 mm (H x W x D)
Weight	2.2 kg (unpacked)
Communications interface	RS232, USB, Bluetooth
Operating conditions	20 - 28°C ; relative humidity: 10 – 85% A dry, clean, level surface. Avoid direct sunlight.
Transport/Storage	0 - 50°C ; relative humidity: 5 – 85%
Power usage	12V DC, 2A

Power Supply

Manufacturer	Globtek
Type	GTM21097-5012
Mains adapter	A separate AC to DC mains adapter with double insulation
Input	90 - 264V AC, 47 - 63 Hz
Output	12V DC, 4.17A
Power usage	Max. 30 VA

SHUTDOWN AND WASTE MANAGEMENT

Before shutting down the MyCare Insite for repair or disposal, ensure that there is no Cuvette left in the Analyser Chamber. To protect the Analyser in shipping the original packaging must be used with the transport lock installed to prevent damage in transport.

For proper waste management according to the Directive 2012/19/EU please contact your Official Saladax Distributor. Used cartridges need to be disposed of with laboratory waste according to the corresponding regulations.

ABOUT THIS USER MANUAL

This user manual will guide you through the installation, operation, and cleaning of your MyCare Insite. The User Manual also explains the operational principles and assists you in troubleshooting any errors or problems. When not used according to the User Manual, the MyCare Insite may not perform properly or be damaged.


For additional information on MyCare Insite functionality and options consult the Technical Guide at [MyCareInsite.com](https://www.MyCareInsite.com).

T +1 (610) 419-6731
F +1 (484) 547-0590
Techsupport@saladax.com

[MyCareInsite.com](https://www.MyCareInsite.com)





© 2024 Saladax Biomedical, Inc.
MyCare™ is a trademark of Saladax Biomedical, Inc.

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MyCare™ is een handelsmerk van Saladax Biomedical, Inc.